



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Antibiotic Zone Reader

Product Image



Description

The Antibiotic Zone Reader provides a rapid accurate method for determining the strength of antibiotic materials by measuring the diameter of an inhibited zone in Petri dish. This Zone Reader measures the diameter of the inhibited zone to 0.02mm with the range of 0 to

35mm works on 230volts A.C. 50 C/S.

-
Antibiotic Zone Reader **Catalog No. 100622**

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.