



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Bottle Top Dispenser, Regular Model

## Product Image



## Description

---

Regular Model :

Wetted parts are PP, Borosilicate Glass (3.3) and FEP. Suitable for routine laboratory application.

### General Product Information:

Description : Regular Model

Piston : PTFE

Cylinder : Borosilicate Glass (3.3)

Volume adjustment Knob : Screw Types

Valve assembly : Glass ball & seat

Valve Housing : Polypropylene

Discharge assembly : PTFE

Delivery tube : FEP

Inlet tube : FEP

Calibration : Individually Calibrated and certified in lab easy calibration by user is also possible.

Accuracy & Reproducibility : In accordance with a ISO 8655 standards

Compatibility : Good for routine lab application

### Regular Model Available in 5 adjustable volume sizes :

**(a)** 0.25-2.5 ml    **(b)** 0.5-5.0 ml    **(c)** 1.0-10ml    **(d)** 2.5-25.0 ml    **(e)** 5.0-50.0ml

- Easy to use, simple construction and smooth, effortless plunger movement for bubble free dispensing.
- Piston : Specially designed PTFE Piston with Silicone O-ring ensures smooth, effortless piston movement and high accuracy.
- Each instrument comes with five adaptors to fit in most of the laboratory reagent bottles.
- Standard threads : 30mm
- Adaptors : 28mm, 32mm, 36mm, 40mm, and 45mm.
- Nozzle cap prevents any unwanted drop on the work space.
- Easy to dis-assemble for cleaning and servicing.
- Glass Barrel is protected by a polypropylene sleeve. It also ensures no glass splitting in case of any accidental breakage.
- Instrument is fully autoclavable at 121 C, 15psi. No Need for disassembling the instrument for autoclaving.
- Screw type volume adjustment knob with large pointer ensures precise, easy and quick volume adjustment.
- Calibration : Specially designed universal calibration tool is provided for convenient

and quick In-lab user re-calibration. This is in compliance with GLO/ISO Norms.

Bottle Top Dispenser, Regular Model

**Catalog No.** 100225

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.