




**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Capillary Tube Apparatus

## Product Image



## Description

---

Ideal for demonstrating the relationship between capillary pressure and the bore diameter of the capillary tube. The apparatus comprises a powder coated metal frame of height 90 mm holding six 150 mm capillary tubes of different bores vertically in a trough. The trough is filled with water and the difference in height of the resulting columns of water is immediately clear.

-

Capillary Tube Apparatus

**Catalog No. 101155**

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.