



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Digital Rotary Microtome

Product Image



Description

It is a modern and reliable sectioning equipment designed for very precise sectioning of tissues upto 1 micron thickness. The feed indicator, located in front of microtome is operated by the latest CAM DRIVE SYSTEM for accurate feed setting from 1-50 microns in steps of one micron. It is graduated in microns, is conveniently located in the front for instant references and its setting is controlled by a knurled knob at the back of the microtomes. The interior mechanism rests on a heavy cast iron base and is covered with a full swing protective cover for easier cleaning and lubrication. The Salient features include independent feed mechanism with automatic safty device, universal knife holder with lateral movements permitting use of the entire knife edge, universal vice type object holder for accurate centring of the specimen.

Particulars

Catalog No.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.