



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Disintegration Test Apparatus (Double Basket)

Product Image



Description

The Apparatus is used for determination of disintegration of time of tablets, a motorized shaft moves in up and down direction at a frequency of 30 cycles per minute, holds 6 glass tubes of standard size with stainless steel wire gauge bottom cylinder which moves up and down with standard dish inside in a liter beaker filled of distilled water which is maintained at 37° C by thermostatically controlled hot plates (without beaker).

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.