



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# E/M by Helical Method Experiment Apparatus

## Product Image

---



## Description

---

**Objective :** To study charge of an electron.

**Features :** Kit Comprises of High Voltage Power Supply Unit with intensity, focus X, Y, deflection & Solenoid current controls.,br> Two metes are provided in Power Supply for Acceleration Voltage control & for solenoid current controls.

One Solenoid Unit Comprises of one Solenoid mounted on Wooden stand & CRT is kept in the solenoid pipe.

<b>Catalog No.</b>	<b>Particulars</b>
JE22110	E/M by Helical Method Experiment Apparatus

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.