



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Expansion of Liquid Apparatus

Product Image



Description

To show difference in expansion of various liquids over the same temperature range. It comprises 5 glass graduated capillary tubes of equal bore and bulbs at the bottom, mounted on a support frame with metal trough for immersion of bulbs in hot water in it. Complete with all accessories including glass parts.

Custom made configuration is available for developing countries or Ministry of Education Tenders

Catalog No.	Particulars
JP10405	Expansion of Liquid Apparatus

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.