



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Manifold Electrophoresis Mechanism

Product Image



Description

It is very ideal to accomplish all types of electrophoresis techniques in very simple way in a single unit. Hence it is a very ideal mechanism for teaching/explain of electrophoresis methods at student level. We can perform the following techniques in this mechanism. 1. Horizontal electrophoresis in gels up to 18x10 cm. 2. Immunoelectrophoresis in 50x50mm or 75x50mm slides. 3. Isoelectric focusing in 18x7cm. gels. 4. Submarine electrophoresis of DNA in 13x7cm. gel. 5. Paper electrophoresis : 3 to 5 strips. 6. Tube gel electrophoresis : 12 tubes. 7. Slab gel electrophoresis : 16x14cm. The manifold electrophoresis mechanism has four compartments :- Horizontal section, Spacer Section, Tube Gel Section & Slab Gel Section. It consists with platinum electrode assembly(2Nos.), Glass Plates for 16x14 cm. gel(1set), Teflon Comb 13 Well, 1.0mm, Spacers 1.0mm(3Nos.), Clamp & Screws(1 set), Glass Tubes 5.0mm I.D. 7.5cm(12Nos.), Rubber Grommets(12Nos.), Gel casting stand,

Immuno Template(1 set), Long needles(2 Nos.), Lid & Power cord

Particulars

Catalog No.

Manifold Electrophoresis
Mechanism

JA20950

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.