



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Micropipette, Regular, Variable Volume

Product Image



Description

Micropipette, Regular is a compact device with separate tip ejector. These micropipettes enjoys excellent reputation with medical researchers and professional people. It comes with superior features to achieve highest efficiency. The unit is user friendly and strain free. Each pipette is calibrated to specifications conforming to DIN 12650 Standards. These pipettes are partially autoclavable.

Features :

- Colour coded for easy identification
- New ergonomic design
- It reduced probability of repetitive strain injury
- Instant volume adjustment
- Separate tip ejector
- Clear digital window on the hand grip

Particulars

Catalog No.

Micropipette, Regular, Variable Volume 0.5-5 ul, Increment 0.1 ul, Accuracy +/-0.1%	JA90500
Micropipette, Regular, Variable Volume 0.5-10 ul, Increment 0.1 ul, Accuracy +/-0.1%	JA90501
Micropipette, Regular, Variable Volume 2-20 ul, Increment 0.1 ul, Accuracy +/-0.1%	JA90502
Micropipette, Regular, Variable Volume 5-40 ul, Increment 0.1 ul, Accuracy +/-0.1%	JA90503
Micropipette, Regular, Variable Volume 5-50 ul, Increment 0.1 ul, Accuracy +/-0.1%	JA90504
Micropipette, Regular, Variable Volume 10-100 ul, Increment 0.1 ul, Accuracy +/-0.1%	JA90505
Micropipette, Regular, Variable Volume 20-200 ul, Increment 0.1 ul, Accuracy +/-0.1%	JA90506
Micropipette, Regular, Variable Volume 40-200 ul, Increment 0.1 ul, Accuracy +/-0.1%	JA90507
Micropipette, Regular, Variable Volume 50-200 ul, Increment 0.1 ul, Accuracy +/-0.1%	JA90508
Micropipette, Regular, Variable Volume 50-500 ul, Increment 0.1 ul, Accuracy +/-0.1%	JA90509
Micropipette, Regular, Variable Volume 100-1000 ul, Increment 0.1 ul, Accuracy +/-0.1%	JA90510
Micropipette, Regular, Variable Volume 200-1000 ul, Increment 0.1 ul, Accuracy +/-0.1%	JA90511

Micropipette, Regular, Variable Volume 1-5 ml, Increment 0.1 ul,
Accuracy +/-0.1%

JA90512

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.