



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Respirometer

Product Image



Description

For comparing rates of absorption or evolution of gases by small organisms, during respiration. The apparatus consists of two stoppered boiling tubes connected by three-way taps to either limb of a graduated manometer. One tube holds the organisms and has a syringe for altering the level of liquid in the manometer during setting up. The other tube acts as a control thermobarometer. CO₂ given off by the respiring organisms is absorbed by potassium hydroxide in the experimental chamber resulting in a change in the level of liquid in the manometer tube.

Respirometer **Catalog No. 100219**

Similar to E8H24890.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.