



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Ria Vial

## Product Image



## Description

---

**Material:** Polypropylene/Polystyrene

These Polypropylene / Polystyrene tubes are generally used in RIA Coagulation & bacteriology. The Polypropylene tubes are clear and can withstand 3000 RCF during centrifugation while the transparent Polystyrene tubes can withstand 1400 RCF. The PP tubes can accept common acids, alkalies & solvents at room temperature while the PS tubes can accept mild bases, weak acids but not organic solvents.

<b>Catalog No.</b>	<b>Particulars</b>
JW06101	12x75mm in PS
JW06102	12x75mm in PP
JW06103	Cap. For 12mm RIA Tube

## Disclaimer

---

The Product details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.