



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Rota Rod Apparatus

Product Image



Description

This Apparatus is latest with digital counter for determining Neuro Pixicity Muscle Tone Balance and Motor Co-ordination in rats and mice consisting of two or three or four compartments of 75mm width each with a rotating rod of 25mm diameter, having speeds 5, 10, 15, 20 and 25 R.P.M. with time interval in each compartment.

-

Rota Rod Apparatus **Catalog No. 100667**

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.