



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Submarine Gel Electrophoresis, Analyst Mechanism

## Product Image



## Description

This is a laboratory scale mechanism. In this mechanism we can run 3 different gels & screen up to 32 samples at a time. This elementary unit consists of detachable platinum electrode assembly(2 nos.), perforated lid for vapors, gel trays 13x13cm., 13x6cm.(2 nos.), 6x6.5cm(4nos.), 15well comb 1.5mm, 8-well comb 1.5mm(2 nos.), 2 well comb 3mm, comb stand connecting cord & instruction manual.

**Particulars**

**Catalog No.**

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.