



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Test Tube Peg Rack

Product Image



Description

Material: Polypropylene

These rack can conveniently hold inverted tubes for drying purpose and minimize collection of airborne contaminants inside the tubes. This rack can also be used to hold and dry electrophoresis and chromatography plates

Catalog No.	Particulars
JW08901	13 mm x 90 Tubes
JW08902	16 mm x 60 Tubes

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.