




Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Thermometer, Laboratory, Superior

Product Image



Description

Glass capillary with bulb type, yellow back with printed scale. Length 300mm with blue plastic ring on top.

Different types available :

Thermometer, Laboratory, Superior, Alcohol in plastic case with blue plastic ring on top

Thermometer, Laboratory, Superior, Mercury in plastic case with blue plastic ring on top

Catalog No.Particulars

JD14680	Thermometer, Laboratory, Superior, Alcohol in plastic case with blue plastic ring on top, -10 to 50 C, 0.5
JD14681	Thermometer, Laboratory, Superior, Alcohol in plastic case with blue plastic ring on top, -10 to 110 C, 1.0
JD14682	Thermometer, Laboratory, Superior, Alcohol in plastic case with blue plastic ring on top, -10 to 110 C, 0.5
JD14683	Thermometer, Laboratory, Superior, Alcohol in plastic case with blue plastic ring on top, -10 to 150 C, 1.0
JD14684	Thermometer, Laboratory, Superior, Alcohol in plastic case with blue plastic ring on top, -50 to + 50 C, 1.0
JD14685	Thermometer, Laboratory, Superior, Alcohol in plastic case with blue plastic ring on top, -50 to + 100 C, 1.0
JD14686	Thermometer, Laboratory, Superior, Alcohol in plastic case with blue plastic ring on top, -100 to + 50 C, 1.0
JD14687	Thermometer, Laboratory, Superior, Alcohol in plastic case with blue plastic ring on top, -100 to + 100 C, 1.0
JD14688	Thermometer, Laboratory, Superior, Alcohol in plastic case with blue plastic ring on top, 0 to 220 F, 1.0
JD14689	Thermometer, Laboratory, Superior, Alcohol in plastic case with blue plastic ring on top, 0 to 400 F, 2.0
JD14690	Thermometer, Laboratory, Superior, Mercury in plastic case with blue plastic ring on top, -10 to 50 C, 0.5
JD14691	Thermometer, Laboratory, Superior, Mercury in plastic case with blue plastic ring on top, -10 to 110 C, 1.0
JD14692	Thermometer, Laboratory, Superior, Mercury in plastic case with blue plastic ring on top, -10 to 110 C, 0.5
JD14693	Thermometer, Laboratory, Superior, Mercury in plastic case with blue plastic ring on top, -10 to 150 C, 1.0
JD14694	Thermometer, Laboratory, Superior, Mercury in plastic case with blue plastic ring on top, -10 to 200 C, 1.0
JD14695	Thermometer, Laboratory, Superior, Mercury in plastic case with blue plastic ring on top, -10 to 250 C, 1.0

- JD14696 Thermometer, Laboratory, Superior, Mercury in plastic case with blue plastic ring on top, -10 to 300 C, 1.0
- JD14697 Thermometer, Laboratory, Superior, Mercury in plastic case with blue plastic ring on top, -10 to 360 C, 2.0

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.