



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Thermometer, Laboratory

Product Image



Description

Glass capillary with bulb type, white/yellow back with printed scale. Length 300mm with ring on top.

Different types available :

Thermometer, Laboratory, Alcohol, White backed, in card board case

Thermometer, Laboratory, Mercury, White backed, in card board case

Thermometer, Laboratory, Alcohol, yellow backed, in card board case

Thermometer, Laboratory, Mercury, yellow backed, in card board case

Thermometer, Laboratory, Alcohol, yellow backed, Deluxe, in plastic case with blue plastic ring on top

Thermometer, Laboratory, Mercury, yellow backed, Deluxe, in plastic case with blue plastic ring on top

Catalog No.	Particulars
JD14620	Thermometer, Laboratory, Alcohol, White backed, in card board case, -10 to 50 C, 0.5
JD14621	Thermometer, Laboratory, Alcohol, White backed, in card board case, -10 to 110 C, 1.0
JD14622	Thermometer, Laboratory, Alcohol, White backed, in card board case, -10 to 110 C, 0.5
JD14623	Thermometer, Laboratory, Alcohol, White backed, in card board case, -10 to 150 C, 1.0
JD14624	Thermometer, Laboratory, Alcohol, White backed, in card board case, -50 to + 50 C, 1.0
JD14625	Thermometer, Laboratory, Alcohol, White backed, in card board case, -50 to + 100 C, 1.0
JD14626	Thermometer, Laboratory, Alcohol, White backed, in card board case, -100 to + 50 C, 1.0
JD14627	Thermometer, Laboratory, Alcohol, White backed, in card board case, -100 to + 100 C, 1.0
JD14628	Thermometer, Laboratory, Alcohol, White backed, in card board case, 0 to 220 F, 1.0
JD14629	Thermometer, Laboratory, Alcohol, White backed, in card board case, 0 to 400 F, 2.0
JD14630	Thermometer, Laboratory, Mercury, White backed, in card board case, -10 to 50 C, 0.5
JD14631	Thermometer, Laboratory, Mercury, White backed, in card board case, -10 to 110 C, 1.0
JD14632	Thermometer, Laboratory, Mercury, White backed, in card board case, -10 to 110 C, 0.5
JD14633	Thermometer, Laboratory, Mercury, White backed, in card board case, -10 to 150 C, 1.0
JD14634	Thermometer, Laboratory, Mercury, White backed, in card board case, -10 to 200 C, 1.0
JD14635	Thermometer, Laboratory, Mercury, White backed, in card board case, -10 to 250 C, 1.0
JD14636	Thermometer, Laboratory, Mercury, White backed, in card board case, -10 to 300 C, 1.0

JD14637	Thermometer, Laboratory, Mercury, White backed, in card board case, -10 to 360 C, 2.0
JD14638	Thermometer, Laboratory, Alcohol, yellow backed, in card board case, -10 to 50 C, 0.5
JD14639	Thermometer, Laboratory, Alcohol, yellow backed, in card board case, -10 to 110 C, 1.0
JD14640	Thermometer, Laboratory, Alcohol, yellow backed, in card board case, -10 to 110 C, 0.5
JD14641	Thermometer, Laboratory, Alcohol, yellow backed, in card board case, -10 to 150 C, 1.0
JD14642	Thermometer, Laboratory, Alcohol, yellow backed, in card board case, -50 to + 50 C, 1.0
JD14643	Thermometer, Laboratory, Alcohol, yellow backed, in card board case, -50 to + 100 C, 1.0
JD14644	Thermometer, Laboratory, Alcohol, yellow backed, in card board case, -100 to + 50 C, 1.0
JD14645	Thermometer, Laboratory, Alcohol, yellow backed, in card board case, -100 to + 100 C, 1.0
JD14646	Alcohol, yellow backed, in card board case, 0 to 220 F, 1.0
JD14647	Thermometer, Laboratory, Alcohol, yellow backed, in card board case, 0 to 400 F, 2.0
JD14648	Thermometer, Laboratory, Mercury, Yellow backed, in card board case, -10 to 50 C, 0.5
JD14649	Thermometer, Laboratory, Mercury, Yellow backed, in card board case, -10 to 110 C, 1.0
JD14650	Thermometer, Laboratory, Mercury, Yellow backed, in card board case, -10 to 110 C, 0.5
JD14651	Thermometer, Laboratory, Mercury, Yellow backed, in card board case, -10 to 150 C, 1.0
JD14652	Thermometer, Laboratory, Mercury, Yellow backed, in card board case, -10 to 200 C, 1.0
JD14653	Thermometer, Laboratory, Mercury, Yellow backed, in card board case, -10 to 250 C, 1.0
JD14654	Thermometer, Laboratory, Mercury, Yellow backed, in card board case, -10 to 300 C, 1.0
JD14655	Thermometer, Laboratory, Mercury, Yellow backed, in card board case, -10 to 360 C, 2.0
JD14656	Thermometer, Laboratory, Alcohol, yellow backed, Deluxe, in plastic case with blue plastic ring on top, -10 to 50 C, 0.5
JD14657	Thermometer, Laboratory, Alcohol, yellow backed, Deluxe, in plastic case with blue plastic ring on top, -10 to 110 C, 1.0
JD14658	Thermometer, Laboratory, Alcohol, yellow backed, Deluxe, in plastic case with blue plastic ring on top, -10 to 110 C, 0.5
JD14659	Thermometer, Laboratory, Alcohol, yellow backed, Deluxe, in plastic case with blue plastic ring on top, -10 to 150 C, 1.0
JD14660	Thermometer, Laboratory, Alcohol, yellow backed, Deluxe, in plastic case with blue plastic ring on top, -50 to + 50 C, 1.0
JD14661	Thermometer, Laboratory, Alcohol, yellow backed, Deluxe, in plastic case with blue plastic ring on top, -50 to + 100 C, 1.0

JD14662	Thermometer, Laboratory, Alcohol, yellow backed, Deluxe, in plastic case with blue plastic ring on top, -100 to + 50 C, 1.0
JD14663	Thermometer, Laboratory, Alcohol, yellow backed, Deluxe, in plastic case with blue plastic ring on top, -100 to + 100 C, 1.0
JD14664	Thermometer, Laboratory, Alcohol, yellow backed, Deluxe, in plastic case with blue plastic ring on top, 0 to 220 F, 1.0
JD14665	Thermometer, Laboratory, Alcohol, yellow backed, Deluxe, in plastic case with blue plastic ring on top, 0 to 400 F, 2.0
JD14666	Thermometer, Laboratory, Mercury, yellow backed, Deluxe, in plastic case with blue plastic ring on top, -10 to 50 C, 0.5
JD14667	Thermometer, Laboratory, Mercury, yellow backed, Deluxe, in plastic case with blue plastic ring on top, -10 to 110 C, 1.0
JD14668	Thermometer, Laboratory, Mercury, yellow backed, Deluxe, in plastic case with blue plastic ring on top, -10 to 110 C, 0.5
JD14669	Thermometer, Laboratory, Mercury, yellow backed, Deluxe, in plastic case with blue plastic ring on top, -10 to 150 C, 1.0
JD14670	Thermometer, Laboratory, Mercury, yellow backed, Deluxe, in plastic case with blue plastic ring on top, -10 to 200 C, 1.0
JD14671	Thermometer, Laboratory, Mercury, yellow backed, Deluxe, in plastic case with blue plastic ring on top, -10 to 250 C, 1.0
JD14672	Thermometer, Laboratory, Mercury, yellow backed, Deluxe, in plastic case with blue plastic ring on top, -10 to 300 C, 1.0
JD14673	Thermometer, Laboratory, Mercury, yellow backed, Deluxe, in plastic case with blue plastic ring on top, -10 to 360 C, 2.0

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.