



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Thin Layer Chromatography Apparatus, Aluminium

## Product Image

---



## Description

---

To Introduce the principles of absorption chromatography using the thin layer technique. This mechanism is having moveable applicator with inbuilt thickness arrangement between 0 to 2 mm (Minimum Div. 0.25 mm). It further consists of the following components. a) Spreader (Applicator) made of ALUMINIUM. b) Perspex base of 114 x 23 cm. to support glass plates. c) Plate rack aluminium, anodised, for ten 20 x 20 cm.

or 20 x 10 cm. plates. d) Developing tank with lid. e) Spotting template made of perspex. f) T.L.C. Plates set of five 20 x 20 cm. and two 20 x 5 cm. or set of ten 20 x 10 cm. and two 20 x 5 cm. g) Glass sprayer with rubber bellow, Cap. 100 ml. h) Micro-pipette. i) Scriber for marking lines made of stainless steel. j) Instruction Manual.

**Particulars****Catalog No.**

Thin Layer Chromatography Apparatus,  
Aluminium

JA20660

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.