



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Tissue Homogeniser

Product Image



Description

It is suitable for low speed homogenisation of tissue and also small volume low speed stirring applications. Low speed operation provides high torque for grinding and homogenising of animal or vegetable tissues. Provided with glass cup with Teflon (P.T.F.E.) pestle having stainless steel shaft, capacity 30 ml and variable speed controlled from 200 to 2400 RPM.

Accessories available at cost :

Teflon Pestle Tissue Homogenisers (Grinder)

It Consists of piston type teflon pestle with serrated tip fitted to a stainless steel rod and a grinding vessel of borosilicate glass for the preparation of tissue homogenates.

Working Capacity :

5 ml.

10 ml.

30 ml.

55 ml.

100 ml.

Teflon Magnetic Rotor

Teflon/P.T.F.E. coated magnetic rotors/needles

Length x Dia.

12 x 8 mm

18 x 8 mm

25 x 8 mm

30 x 8 mm

35 x 8 mm

50 X 9 mm

Teflon Policeman 300 mm

Tissue Homogeniser

Catalog No. 100235

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.