



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

VDRL Shaker

Product Image



Description

It is very useful for V.D.R.L. grouping & typing as well as other agglutination tests. It is also used for mixing of solution in small flask, beakers, bottles or similar containers. It has a timer incorporated in the electrical circuit which is being cut off automatically as soon as the setting is over. Platform 30x30 cm. Also available with digital timer.

Particulars	Catalog No.
V.D.R.L. Shaker	JA 23830
V.D.R.L. Shaker with Digital RPM & Timer	JA 23831

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.