



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Tissue Culture Rack

## Product Image

---



## Description

---

**Construction:** Tissue racks are made of tubular mild steel pipes duly painted for longer life/ made of aluminium.

**Racks:** Each rack of size 48"x12" (depth) is covered with thick unbreakable plexi glass/sunmica paste duro board, which are terminate and insect proof fire and weather resistant. Distance between each rack is 15".

**Lighting:** Four fluorescent tubes are properly illuminated by diffused, glare free

fluorescent lights are built in upper portion in worktable. UV Germicidal Light are fitted inside the rack/ work area. All fluorescent tubes chocks and ultra violet germicidal tube chocks are fitted outside of rack. If desired to avoid heat dissipation inside the rack.

**Timer:** 0-24 Hours Digital Timers can be provided in tissue culture rack.

**Control Panels:** All control panels are fitted tops of tissue culture racks including On/Off switch of fluorescent tubes and UV tubes. All control panels are fitted separately in each rack with indicator light.

**Easy Mobility:** The whole rack is provided with four castor wheels at the bottom for easy mobility to any place.

Unit of four shelves with 16 Fluorescent tubes 4 U.V. lights.

Unit of five shelves with 20 Fluorescent tubes 5 U.V. lights.

Unit of six shelves with 24 Fluorescent tubes 6 U.V. lights.

Tissue Culture Rack

**Catalog No.** 100398

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.